

DSLR camera system...



Features include:

- Live view on PC, &/or LCD/Plasma TV
- High-resolution still images (12.1MP to 18.1MP)
**depending on camera model
- HD Video recording (720P or 1080P)
- Compatible with most DSLR cameras (Canon, Nikon, Olympus, Fuji, Minolta, Sony (A-series, and NEX 3/5/7))
- Wireless download capable using Eye-fi card (www.eye.fi)
- Imaging software also available...

Easy Setup



- Remove Binocular head
- Mount beamsplitter above, below or to the side of the binoculars.
(beamsplitter available in single port vertical mount or single/dual port side mount, depending on model of slit lamp)
- Replace Binocular head
- Mount DSLR Camera adaptor
- Mount DSLR camera to camera adaptor

Compatible with the following Slit Lamps and Surgical microscopes...

Appasamy A1A11-SS; Burton 860; Keeler SL – 33, 40, 40HT; Haag Streit 900BQ; Marco G2, G4 & G5; Nidek SL 1600/1800; Reichert Xcel 250,255, 400 & 700; Shin Nippon SL 203, 0450X, D2, D4, D7 & D8Z; ALL Topcon slit lamps and OMS microscopes (**except for SL 3 series); Veatch ReSeevit 3Z, 5V & 5Z; Woodlyn HR-3, HR-5, SMM-3, SMM-5; ALL Leica surgical microscopes; ALL MOLLER-WEDEL surgical microscope, ALL Zeiss slit lamps and Microscopes.



2246 Camino Ramon • San Ramon • CA 94583 • USA
Tel: 925-355-0750 • 800-322-7373 • Fax: 925-355-0777
e-mail: info@ttimedical.com • web site: www.ttimedical.com

